


Technical Sheet

# INTEGRALE<sup>®</sup> & GENERIC<sup>®</sup>

Range of Straight Femoral Stems  
Cemented and Cementless



AMPLITUDE<sup>®</sup>


# INTEGRALE® & GENERIC®

## INTEGRALE® & GENERIC®: COMMON FEATURES

- Straight femoral stems product line
- Ovoid cross-section for an optimized metaphyseal filling
- Neck-shaft angle: 137°
- 7 sizes : homothetic increase in sizes
- Small Morse Taper 10/12 (5°42'30'')
- Mirror-polished thin neck
- Same intra-medullary features for both stems
- Same stem length for both stems
- Rounded edges
- Engraving on proximal part of Morse Taper
- Oblong hole for impaction and threaded hole for extraction, positioned at the top of the stem



## INTEGRALE® : SPECIFIC FEATURES




- Perfect visualization of the impaction limit : HA coating
- Longitudinal grooves closed on proximal side
- HA coating 80µm
- Material: Titanium Alloy(Ti6Al4V) anodized

## GENERIC® : SPECIFIC FEATURES



- Highly-polished surface on the intra-medullary part of the stem
- Material: Stainless Steel (M30NW)

		INTEGRALE® GENERIC®	LENGTH L
Size 1	O	41.2	111
	H	34.3	
Size 2	O	41.7	116
	H	34.3	
Size 3	O	42.2	122
	H	34.3	
Size 4	O	42.7	127
	H	34.3	
Size 5	O	43.2	133
	H	34.3	
Size 6	O	43.7	138
	H	34.3	
Size 7	O	44.2	144
	H	34.3	



Nearly 20 years  
of clinical  
experience

\*Distance between diaphysis and head centre with medium neck



The INTEGRALE® system is available on AMPLIVISION®:  
Assists the surgeon for extra-medullary settings.

**AMPLIVISION®**

# References



## INTEGRALE®

INTEGRALE® femoral stem - cementless - HA coated - Size <b>1</b>	1-01001 <b>01</b>
INTEGRALE® femoral stem - cementless - HA coated - Size <b>2</b>	1-01001 <b>02</b>
INTEGRALE® femoral stem - cementless - HA coated - Size <b>3</b>	1-01001 <b>03</b>
INTEGRALE® femoral stem - cementless - HA coated - Size <b>4</b>	1-01001 <b>04</b>
INTEGRALE® femoral stem - cementless - HA coated - Size <b>5</b>	1-01001 <b>05</b>
INTEGRALE® femoral stem - cementless - HA coated - Size <b>6</b>	1-01001 <b>06</b>
INTEGRALE® femoral stem - cementless - HA coated - Size <b>7</b>	1-01001 <b>07</b>

## GENERIC®

GENERIC® femoral stem - Cemented - size <b>1</b>	1-01002 <b>01</b>
GENERIC® femoral stem - Cemented - size <b>2</b>	1-01002 <b>02</b>
GENERIC® femoral stem - Cemented - size <b>3</b>	1-01002 <b>03</b>
GENERIC® femoral stem - Cemented - size <b>4</b>	1-01002 <b>04</b>
GENERIC® femoral stem - Cemented - size <b>5</b>	1-01002 <b>05</b>
GENERIC® femoral stem - Cemented - size <b>6</b>	1-01002 <b>06</b>
GENERIC® femoral stem - Cemented - size <b>7</b>	1-01002 <b>07</b>

## Femoral heads

	Short neck	Medium neck	Long neck	Extra-long neck
Metallic femoral head M30NW - Taper 10/12 - Ø <b>22.2 mm</b>	1-01011 <b>01</b>	1-01011 <b>02</b>	1-01011 <b>03</b>	-
Metallic femoral head M30NW - Taper 10/12 - Ø <b>28 mm</b>	1-01011 <b>11</b>	1-01011 <b>12</b>	1-01011 <b>13</b>	1-01011 <b>14</b>
Metallic femoral head M30NW - Taper 10/12 - Ø <b>32 mm*</b>	1-01011 <b>21*</b>	1-01011 <b>22*</b>	1-01011 <b>23*</b>	-
BIOLOX® Delta ceramic femoral head - Taper 10/12 - Ø <b>28 mm</b>	1-01911 <b>01</b>	1-01911 <b>02</b>	1-01911 <b>03</b>	-
BIOLOX® Delta ceramic femoral head - Taper 10/12 - Ø <b>32 mm</b>	1-01911 <b>11</b>	1-01911 <b>12</b>	1-01911 <b>13</b>	-
BIOLOX® Delta ceramic femoral head - Taper 10/12 - Ø <b>36 mm</b>	1-01911 <b>21</b>	1-01911 <b>22</b>	1-01911 <b>23</b>	-

\*optional, available on request only

**Service Clients-France :**  
 Porte du Grand Lyon  
 01700 Neyron, France  
 Tel. : +33 (0)4 37 85 19 19  
 Fax : +33 (0)4 37 85 19 18  
 E-mail : [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

**Customer Service-Export :**  
 11, cours Jacques Offenbach, ZA Mozart 2  
 26000 Valence, France  
 Tel. : +33 (0)4 75 41 87 41  
 Fax : +33 (0)4 75 41 87 42  
 E-mail : [amplitude@amplitude-ortho.com](mailto:amplitude@amplitude-ortho.com)

Internet : [www.amplitude-ortho.com](http://www.amplitude-ortho.com)