

SCORE® AS (Allergy Solution)



- SCORE® AS Total knee prosthesis coated with Titanium Nitride (TiN) ceramic to minimise allergic reactions
- Coating acts as a barrier to reduce the release of allergenic ions
- Coating thickness: 4 µm
- Cemented components:
 - Femoral component
 - Tibial base plate
 - Tibial extension stems
- Highly polished articular surfaces:
 - Femoral component / tibial insert
 - Tibial Insert / tibial base plate
- Available only upon request

SCORE® Allergy Solution femoral component PCL sacrificing cemented Size 1 to Size 7 RIGHT	1-0206101 to 1-0206107
SCORE® Allergy Solution femoral component PCL sacrificing cemented Size 1 to Size 7 LEFT	1-0206111 to 1-0206117
SCORE® Allergy Solution tibial base plate for mobile bearing insert cemented Size 1 to Size 7	1-0206201 to 1-0206207
SCORE® mobile bearing insert - Size 1 to 7 thickness 10 mm	1-0200611 to 1-0200617
SCORE® mobile bearing insert - Size 1 to 7 thickness 12 mm	1-0200621 to 1-0200627
SCORE® mobile bearing insert - Size 1 to 7 thickness 14 mm	1-0200631 to 1-0200637
SCORE® mobile bearing insert - Size 1 to 7 thickness 16 mm	1-0200641 to 1-0200647
SCORE® mobile bearing insert - Size 1 to 7 thickness 20 mm	1-0200651 to 1-0200657

Available upon request

Allergy Solution extension stem Ø 10 length 75 mm	1-0207310
Allergy Solution extension stem Ø 10 length 100 mm	1-0207311
Allergy Solution extension stem Ø 12 length 75 mm	1-0207320
Allergy Solution tibial half block Size 1/2 thickness 5 mm	1-0207511
Allergy Solution tibial half block Size 3/4/5 thickness 5 mm	1-0207513
Allergy Solution tibial half block Size 6/7 thickness 5 mm	1-0207516