



TECHNICAL SHEET SURGICAL CEMENT

 **amplifix 1**

 **amplifix 1G**
With gentamicine

 **amplifix 3G**
With gentamicine

AMPLIFIX 1


Name	Product No.
Amplifix 1 - Surgical cement standard viscosity	1-0400101

Indication : AMPLIFIX 1 bone cement is indicated for fixation of prostheses to living bone in orthopaedic musculoskeletal surgical procedures for rheumatoid arthritis, osteoarthritis, traumatic arthritis, osteoporosis, avascular necrosis, collagen disease, severe joint destruction secondary to trauma or other conditions, and revision of previous arthroplasty.

Medical Device class : IIb

AMPLIFIX 1 and 3 - Gentamicin


Name	Product No.
Amplifix 1 G- Surgical cement high viscosity	1-0400201
Amplifix 3 G - Surgical ciment low viscosity	1-0400203

Indication : AMPLIFIX bone cements with gentamicin are intended for fixation of prosthetic components into bone medullar cavity in arthroplasty procedures.

Medical Device class : III

Packaging

MONOMER : sterile liquid ampoule in a blister

POLYMER : sterile bag powder (40 grams for the version without gentamicin and 41,6 grams for the version with gentamicin)

Double sterile packaging

STORAGE : maximum temperature of 25 °C

Legal Manufacturer : TEKNIMED SAS, 8, rue du Corps Franc-Pommiès, 65500 VIC en BIGORRE, France

Distributed by AMPLITUDE

Customer Service – France :

Porte du Grand Lyon,
01700 Neyron – France
Tél. : +33 (0)4 37 85 19 19
Fax : +33 (0)4 37 85 19 18
E-mail : amplitude@amplitude-ortho.com

Customer Service – Export :

11, cours Jacques Offenbach,
ZA Mozart 2,
26000 Valence – France
Tél. : +33 (0)4 75 41 87 41
Fax : +33 (0)4 75 41 87 42