


INTEGRALE[®] Revision cementless
GENERIC[®] Revision cemented

Revision femoral stems



AMPLITUDE[®]

INTEGRALE® & GENERIC® Revision


INTEGRALE® & GENERIC® REVISION: COMMON FEATURES

- Straight femoral stems product line
- Ovoid cross-section for an optimized metaphyseal filling
- Neck length: + 3,5 mm compared to primary stems of the INTEGRALE® range
- Neck-shaft angle: 135°
- 7 sizes : homothetic increase in sizes
- Small Morse Taper 10/12 (5°42'30'')
- Mirror-polished thin neck
- Same intra-medullary features for both stems
- Same stem length for both stems
- Rounded edges
- Engraving on proximal part of Morse Taper
- Oblong hole for impaction and threaded hole for extraction, positioned at the top of the stem



INTEGRALE® REVISION: SPECIFIC FEATURES

- Perfect visualization of the impaction limit : HA coating
- Longitudinal grooves closed on proximal side
- HA coating 80µm
- Material: Titanium Alloy (Ti6Al4V) anodized



GENERIC® REVISION: SPECIFIC FEATURES

- Highly-polished surface on the intra-medullary part of the stem
- Material: Stainless Steel (M30NW)




		INTEGRALE® GENERIC®	LENGTH L
Size 1	O	41.2	164
	H	37.5	
Size 2	O	41.7	169
	H	37.5	
Size 3	O	42.2	174
	H	37.5	
Size 4	O	42.7	179
	H	37.5	
Size 5	O	43.2	184
	H	37.5	
Size 6	O	43.7	189
	H	37.5	
Size 7	O	44.2	194
	H	37.5	



Nearly 20 years
of clinical
experience

*Distance between diaphysis and head centre with medium neck



References

CE₂₇₉₇

INTEGRALE® REVISION

INTEGRALE® Revision femoral stem - HA coated - cementless - Size 1	1-01003 01
INTEGRALE® Revision femoral stem - HA coated - cementless - Size 2	1-01003 02
INTEGRALE® Revision femoral stem - HA coated - cementless - Size 3	1-01003 03
INTEGRALE® Revision femoral stem - HA coated - cementless - Size 4	1-01003 04
INTEGRALE® Revision femoral stem - HA coated - cementless - Size 5	1-01003 05
INTEGRALE® Revision femoral stem - HA coated - cementless - Size 6	1-01003 06
INTEGRALE® Revision femoral stem - HA coated - cementless - Size 7	1-01003 07

GENERIC® REVISION

GENERIC® Revision femoral stem - Cemented - Size 1	1-01004 01
GENERIC® Revision femoral stem - Cemented - Size 2	1-01004 02
GENERIC® Revision femoral stem - Cemented - Size 3	1-01004 03
GENERIC® Revision femoral stem - Cemented - Size 4	1-01004 04
GENERIC® Revision femoral stem - Cemented - Size 5	1-01004 05
GENERIC® Revision femoral stem - Cemented - Size 6	1-01004 06
GENERIC® Revision femoral stem - Cemented - Size 7	1-01004 07

Femoral heads

	Short neck	Medium neck	Long neck	Extra-long neck
Metallic femoral head M30NW - Taper 10/12 - Ø 22.2 mm	1-01011 01	1-01011 02	1-01011 03	-
Metallic femoral head M30NW - Taper 10/12 - Ø 28 mm	1-01011 11	1-01011 12	1-01011 13	1-01011 14
Metallic femoral head M30NW - Taper 10/12 - Ø 32 mm*	1-01011 21*	1-01011 22*	1-01011 23*	-
BILOX® Delta ceramic femoral head - Taper 10/12 - Ø 28 mm	1-01911 01	1-01911 02	1-01911 03	-
BILOX® Delta ceramic femoral head - Taper 10/12 - Ø 32 mm	1-01911 11	1-01911 12	1-01911 13	-
BILOX® Delta ceramic femoral head - Taper 10/12 - Ø 36 mm	1-01911 21	1-01911 22	1-01911 23	-

Service Clients-France :
Porte du Grand Lyon
01700 Neyron, France
Tel. : +33 (0)4 37 85 19 19
Fax : +33 (0)4 37 85 19 18
E-mail : amplitude@amplitude-ortho.com

Customer Service-Export :
11, cours Jacques Offenbach, ZA Mozart 2
26000 Valence, France
Tel. : +33 (0)4 75 41 87 41
Fax : +33 (0)4 75 41 87 42
E-mail : amplitude@amplitude-ortho.com

Internet : www.amplitude-ortho.com

AMPLITUDE®