



TECHNICAL SHEET



Range of «taper slip» cemented femoral stems

Up to 4 offsets depending on size

Based on a design with 50 years of clinical evidence

General design



13 stems available in the range, up to 4 offsets per size, up to 5 sizes per offset

Material: Stem = Stainless steel (M30NW)

Distal centralizer = PMMA

Important Features

Instruments



Oversized broaches for cement mantle thickness



Ergonomic gun for stem insertion:
prehension and release of the stem
in one gesture

Dimensions and size chart

Size 1: Offset 35,5 - Offset 37,5 - Offset 44 - Offset 50 mm



Size	Length L (mm)	Offset O (mm)				Neck Length C (mm)	Neck-shaft Angle
0		-	37,5	44	-		
1	130 to	35,5	37,5	44	50		
2	134	-	37,5	44	50	32,3	125°
3	(107 for Size 1 offset 35,5)	-	37,5	44	50		
4		-	-	44	-		

E² FEMORAL STEM	35,5 mm Offset	37,5 mm Offset	44 mm Offset	50 mm Offset
E ² Femoral Stem - Cemented - Taper 12/14 - Size 0	-	111-042-005	111-042-009	-
E ² Femoral Stem - Cemented - Taper 12/14 - Size 1	111-042-002	111-042-006	111-042-010	111-042-014
E ² Femoral Stem - Cemented - Taper 12/14 - Size 2	-	111-042-007	111-042-011	111-042-015
E ² Femoral Stem - Cemented - Taper 12/14 - Size 3	-	111-042-008	111-042-012	111-042-016
E ² Femoral Stem - Cemented - Taper 12/14 - Size 4	-	-	111-042-013	-

FEMORAL HEADS	SHORT NECK	MEDIUM NECK	LONG NECK	EXTRA-LONG NECK
BIOLOX® Delta ceramic femoral head - Taper 12/14 — Ø28 mm	1-0190911	1-0190912	1-0190913	-
BIOLOX® Delta ceramic femoral head - Taper 12/14 — Ø32 mm	1-0190921	1-0190922	1-0190923	1-0190924
BIOLOX® Delta ceramic femoral head - Taper 12/14 — Ø36 mm	1-0190931	1-0190932	1-0190933	1-0190934
CoCr Metal Femoral Head - Taper 12/14 — Ø22,2 mm	111-022-001	111-022-002	111-022-003	-
CoCr Metal Femoral Head - Taper 12/14 — Ø28 mm	111-028-001	111-028-002	111-028-003	111-028-004
CoCr Metal Femoral Head - Taper 12/14 — Ø32 mm	111-032-001	111-032-002	111-032-003	111-032-004
CoCr Metal Femoral Head - Taper 12/14 — Ø36 mm	111-036-001	111-036-002	111-036-003	111-036-004

Service Clients – France :

 Porte du Grand Lyon,
 01700 Neyron – France

 Tél. : **+33 (0)4 37 85 19 19**

Fax : +33 (0)4 37 85 19 18

 E-mail : amplitude@amplitude-ortho.com
Customer Service – Export :

 11 cours Jacques Offenbach,
 ZA Mozart 2,

26000 Valence – France

 Tél. : **+33 (0)4 75 41 87 41**

Fax : +33 (0)4 75 41 87 42

www.amplitude-ortho.com